


5 Easy Steps to X-Ray Identification



1 Cut tape to desired length or use our convenient pre-cut, 2 and 3-inch labels.



2 Write or type patient data on label.



3 Remove label from the easy release crack-and-peel backing paper.



4 Adhere label to your personalized density filter.




5 Place filter on outside of film cassette. Expose the film. Patient and practice data becomes a permanent, clear, high contrast image on your x-rays.

MED I.D.[®]

Medical I.D. Systems Inc.
3954 44th Street S.E.
Grand Rapids, MI 49512

BULK RATE
U.S. POSTAGE
PAID
PERMIT NO.
657
GRAND RAPIDS, MI




MED I.D.[®]

X-Ray Label Tape

Saves Time, Saves Money

Eliminate X-Ray Misidentification Risks. . .
SEE INSIDE!



“The clarity and readability of MED I.D.[®] x-ray identification tape is excellent.”

Harriet White,
PMAC

Eliminate X-Ray Misidentification Risks

. . . SEE INSIDE!

